

## Delta® Ultra-Light 1000 Full-Electric Low Bed

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- The Delta 1000 is truly universal because the bed frame can be used with Drive's and most other manufacturers' bed ends, old and new
- The transition box mounted on the foot section is key to interchangeability (Figure A)
- The transition box allows you to change the rotation of the high/low shaft
- Headboard is taller than footboard to better fit home decor (Figure B)
- The wood grain panels are break and scratch resistant
- Once installed, the high/low shaft will stay with the bed to eliminate loss or forgetting on future installations (Figure C)
- Labels and color-coded springs assure accurate installation of side rails for patient safety (Figure D)
- Unique motor is completely self-contained to reduce weight and noise
- 9 ½" to 23 ½" deck height
- Easy to set up
- Top mounted Hi/Lo crank accessibility
- Height-adjustable motor comes pre-installed on the foot section
- Caster easily moves from standard to low bed position
- The motor assembly can be installed and removed with patient in bed
- Bed, side rails and mattress safety-tested by SGS, a nationally recognized testing facility Satisfies all applicable FDA entrapment guidelines
- The motor is UL approved
- In the event of a power failure, a 9 volt battery, located in the hand control, can lower the head and foot sections nine times
- New and improved hand pendant has large, easy-to-use controls (Figure E)
- Bed ships in two cartons
- Hand crank located on the top of the footboard prevents bending for the caregiver when in the low position (Figure F)
- Channel frame construction provides superior strength and reduced weight
- Zinc-coated spring deck
- Extension kit available (model #15030EXTKIT) extends length from 80" to 84"
- Patent Pending
- Weight Capacity 450 lbs

PRODUCTS	Description
15235	LowBed,UltraLightFull.Elec.1/c
15235BV-FR	LowBed,UltraLightFull.Elec.1/c
15235BV-HR	LowBed,UltraLightFull.Elec.1/c
15235BV-PKG	LowBed,UltraLightFull.Elec.Pkg
15235BV-PKG-1	LowBed,UltraLightFull.Elec.Pkg
15235BV-PKG-1-T	LowBed,UltraLightFull.Elec.Pkg
15235BV-PKG-T	LowBed,UltraLightFull.Elec.Pkg

SPECIFICATIONS	Description
Warranty	(Original Purchaser) 5 Years on Motor and Hand Control, Limited Lifetime on Welds and Frame, 1 Year on All Other Parts and Components
Weight Capacity	450 lbs Maximum weight limit includes weight of occupant
Carton Shipping Weight	191 lbs.
Weight	175 lbs.
ACCESSORIES	Description
15030HBL	Hospital Bed Linen Kit,1 EA,10/CS, CS
15030HBC	Bedding in a Box,1 EA,5/CS, CS
PARTS	Description
15044CSET	Set of 4 Casters,Homecare Bed, 1 set/ea, EA
15235	LowBed,UltraLightFull.Elec.1/c, CS
15547F-FP-1	Foot Board Panel, 15235, Delta Bed, 1/ea, EA
15547H-HP-1	HeadBoardPanel,15235/15547H,DeltaBed,1ea, EA
15044CL	Caster,w/Lockfor15544,15547,1/ea, EA
15044CNL	Caster,NoLockfor15544,15547,1/ea, EA
15541	Head Spring,15030/15033,Delta Bed,1ea, EA
15030-FRLINK	FrameLinkw/DeckBolts,15030,15005, 1/ea, EA
15544HGB	Head Board Gear Box,15544,1ea, EA
15544FGB	Foot Board Gear Box,15544,1ea, EA
1503C	Locking Collar Universal Bed Shaft,1ea, EA
15544GBMB	MountingBracketForTransBox15030,1ea, EA
15033MHLNC-N	HiLoMotorConnector15033MHLN-MGBA,1ea, EA
15030MC	Cap Only for 15030 Motor,1/ea, EA
15030-MS	Motor Strap,Ultra Lite Bed,2 per PKG, PKG
15030PN-1035	HeadActuatorBarDeltaFlatTubing,1ea, EA
15030PN-1036	FootActuatorBarDeltaFlatTubing,1ea, EA
15030-7RMM	Motor Mount,7 R Serial,1/ea, EA
15030P-1006	Covers,UltraLt.Bed Okin,1/bx, BX
15538	Head/Foot Board, 15235N, 1/ea, EA