

22" Bariatric Aluminum Transport Chair

HCPCS: E1039



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E




FIGURE F



FIGURE G

- Weighs only 33 lbs
- Supports individuals up to 450 lbs
- Attachable push bar with padded foam grip
- 12" rear flat-free wheels
- Back folds down for convenient storage
- Padded armrests provide additional comfort
- Heavy-duty, nylon reinforced upholstery
- Comes with hand wheel lock levers and companion activated wheel locks
- Comes with seat belt
- Comes with swing-away footrests and anti-tippers

PRODUCTS

Description

ATC22-R	Aluminum HD 22" Red Transport Chair,1/ea
---------	--

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Product Dimensions	41" (L) w/ riggings x 30" (W) x 37.5" (H)
Seat Dimensions	22" (W) x 18" (D)
Weight Capacity	450 lb
Seat To Floor Height	21"
Carton Shipping Weight	39 lb
Overall Width	30" (Open) 12" (Closed)
Weight	33 lb
Back Height	16"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]