



PROUDLY MADE  
IN THE U S A

LIMITED  
LIFETIME  
WARRANTY



SAFETY | INTEGRITY



SECURITY | QUALITY

**EMPOWERING  
INDEPENDENCE™**



## ASSURANCE™ with Clean-Shield™



### 3" RAISED PREMIUM TOILET SEAT + OPTIONS

- Heavy-duty, Oversized, Raised Seat for Comfort & Support
- Seat Capacity: 1,000 lbs.
- Clean-Shield™ Keeps Toilet Cleaner
- Stay-Tite™ Commercial Fastening System™ Never Loosens. Guaranteed.
- Larger Seat Opening
- Non-corroding Stainless Steel Posts & Pintles
- Four Super-Grip™ Bumpers
- Self-sustaining Check Hinges Keep Seat Upright for Standing Use or Cleaning
- Easy-to-Clean Plastic with Residential Style Cover



#### ASSURANCE™ (seat only)


Elongated | E05310TSS  
Round | R05310TSS



#### + ASSURANCE™ SUPPORT ARMS

A05313 GRY  
(arms only)

- Ergonomic Design Better Supports Sitting and Standing
- 21" between Support Arms
- Arm Capacity: 350 lbs. each



#### + ASSURANCE™ BIDET ATTACHMENT\*

A06303TWA  
(attachment only)

- Non-electric, Simple Twist Knob Controls Wash During Use
- Adjustable T-valve for Customizable and Consistent Water Pressure
- Water Spray Angle Accommodates Higher Seat Height



\* Non-domestic




## ASSIST™



### STANDARD HEIGHT PREMIUM TOILET SEAT + SUPPORT ARMS

- Replaces Need for Wall- or Toilet-mounted Rails
- Seat Capacity: 500 lbs. | Arm Capacity: 250 lbs. each
- Stay-Tite™ Commercial Fastening System™ Never Loosens. Guaranteed.
- Larger 8-1/4 x 12-1/8" Seat Opening
- Seat Lifts for Standing Use & Cleaning
- 12 Stabilizing Seat Ribs
- Easy-to-Clean Plastic with Residential Style Cover
- 21-3/4" between Support Arms (Included)



Universal-Fit | A0AS100  
Fits Both Round & Elongated Toilets




## RISE™



### 4.5" BASIC RAISED TOILET SEAT + SUPPORT ARMS

- Taller Height for Additional Support
- Seat Capacity: 300 lbs. | Arm Capacity: 200 lbs. each
- Stay-Tite™ Seat Fastening System™ Never Loosens. Guaranteed.
- Larger 8-1/4 x 11-3/4" Seat Opening
- Easy-Clean® - Seat Removes for Easy Cleaning
- Easy-to-Clean Plastic
- Coverless | Fixed Mount (Does not lift for standing use)
- 18-9/16" between Support Arms (Included)



Universal-Fit | A04505T  
Fits Both Round & Elongated Toilets  
Ideal for Temporary Use





## STEADFAST™

## 3" RAISED OPEN-FRONT HEAVY-DUTY TOILET SEAT



Elongated | E82350TC

- Raised Seat for Comfort & Support
- Stay-Tite™ Commercial Fastening System™ Never Loosens. Guaranteed.
- Non-corroding Stainless Steel Posts & Pintles
- Heavy-duty Commercial Grade Easy-to-Clean Plastic
- Open-Front Provides Easier Access
- Four Super-Grip™ Bumpers
- Self-sustaining Check Hinges Keep Seat Upright for Standing Use & Cleaning
- Seat Capacity: 500 lbs.



# STAY·TITE™

COMMERCIAL FASTENING SYSTEM™

## NEVER LOOSENS. GUARANTEED.

**STRONG**

Superior, reinforced,  
& non-corroding

**FAST**

Easy & precise  
installation

**SECURE**

Never re-tighten  
seat again

HARDWARE  
NEVER  
LOOSENS



Nut snaps off at  
perfect torque, so  
it can't be over or  
under-tightened.



Family-owned since 1901, Bemis Manufacturing Company has been improving customers' lives by engineering seats with trust-worthy features and quality craftsmanship made in the USA. Our residential style seats have features that simplify life, providing stability, comfort and cleanliness. We take great pride in ensuring both value and quality in our innovative products for everyday life.



Bemis Manufacturing Co.  
1.888.722.6488 | ToiletSeats.com

©2023 0B7012562