



BESTLIFT SERIES

FULL BODY PATIENT LIFT

SPECIALTY FULL BODY PATIENT LIFT

PL350CT [Car Transfer](#)
PL400EF [Foldable/Portable](#)
PL 400H/HE [Convertible](#)

CARE FULL BODY PATIENT LIFT

PL 182
PL 228

PROCARE FULL BODY PATIENT LIFT

PL 400
PL 500
PL 600

BESTLIFT®

PL350CT | FULL BODY PATIENT LIFT

The BestLift® PL350CT is the perfect solution for car transfers and your everyday transfer needs. With its unique spreader bar, this BestLift is able to position patients from sitting to a lying down position with minimal patient movement. It also transfer them from the car, wheelchair, floor or bed and adds versatility by folding up for easy transportation.

Car Transfer Specialist

- **350lb** capacity
- Easy car transfer
- Fully electric lift with optimal spare battery
- Foldable and portable for easy transport
- Compatible with most plastic clip slings
- Minimal swinging of patient



2920 PACIFIC DR., NORCROSS, GA 30071
WWW.BESTCARELLC.COM
678-679-6690



| PL350CT - Car Transfer Lift



PL350CT BEING USED TO
PICK UP A PATIENT



SIDE VIEW OF THE PL350CT
INCLINED POSITION



SIDE VIEW OF THE PL350CT
LYING POSITION


Easy Transfer in-and-out cars




| DIMENSIONS

Key dimensions of PL350CT	
maximum capacity	350Lb / 158.8Kg
base legs length	41.0 in
base closed external width	25.6 in
base opened internal width	32 in
base height	4.7 in
mast height	47.3 in
min lifting height	21.5 in
max lifting height	63.5 in
netweight	94Lb / 43Kg
shipping weight	108Lb / 49Kg

SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE
AND ARE SUBJECT TO CHANGE WITHOUT NOTICE



TOP VIEW OF PL350CT



SIDE VIEW OF PL350CT

BESTLIFT®

NEW!

PL400EF | FULL BODY PATIENT LIFT

The Bestlift® PL400EF represents a new standard in full body lifts for individuals who travel or are seeking storage efficiency in their home. PL400EF assembles (or folds) with 3 easy steps allowing for full operation or preparing for storage/transport.




PUSH HANDLES



- Portable and foldable for easy storage and transportation
- Light weight and easy to move
- **400lb** capacity
- Optional scale available
- Optional spare battery and wall charger available
- Easy to service and disassemble


| PL400EF




| DIMENSIONS

Dimensions & Specifications		
	Safe Working Load	400 lb (182 kg)
B	Base Length	41.0" (1,041mm)
C	Base External Width – Closed	25.6" (650mm)
D	Base Internal Width – Opened	32.0" (813mm)
E	Mast Height	48.3" (1,227mm)
F	Front Caster Height	4.7" (119mm)
G	Lifting Height – Lowest	18.5" (465mm)
H	Lifting Height – Highest	67.2" (1707mm)
I	Base Opened Turning Diameter	55" (1,384.9mm)
H	Base Closed Turning Diameter	51" (1,306.5mm)
	Unloaded Lift Weight	89.7 lb (40.8 kg)

* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE
AND ARE SUBJECT TO CHANGE WITHOUT NOTICE





TOP VIEW OF PL400EF



SIDE VIEW OF PL400EF





2920 PACIFIC DR., NORCROSS, GA 30071
WWW.BESTCARELLC.COM
678-679-6690



BESTLIFT™

PL400H/HE | CONVERTIBLE

The innovative BestLift™ PL400H/HE home care lifts are designed to allow manual or electric use. The hydraulic unit may be upgraded to electric anytime with an easy-to-attach electric kit. Built to last and yet extremely affordable. Ask your local dealer about the Medicare options available for this lift.



- 6 point spreader bar
- **400lb** capacity
- Electric
- Medicare options available
- No tool required for assembly

BestLift PL400HE Displayed



BestLift PL400H Displayed

- 6 point spreader bar
- **400lb** capacity
- Hydraulic
- Medicare options available

| UPGRADE ELECTRONIC KIT



| SPECIAL FEATURE

- *EASY TO DISASSEMBLE
FOR TRANSPORTING
- *NO TOOLS NEEDED



| DIMENSIONS

Key dimensions of PL400H/HE

maximum capacity	400Lb / 182Kg
base legs length	40"
base closed external width	26"
based open external width	33"
mast height	49"
base height	4.8"
min. lifting height	27"
max. lifting height	68"
net weight	82Lb / 37Kg
shipping weight	93Lb / 42Kg

* with spreader bar

SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE
AND ARE SUBJECT TO CHANGE WITHOUT NOTICE



BESTLIFT™

NEW!

PL182 | FULL BODY PATIENT LIFT

The BestLift™ PL182 is the perfect blend of form and function for patient transfers. With an overall base length of 44.9", maximum width of 42.9", and base height of 4.5", it is the perfect solution for transferring bariatric patients in small rooms with low beds.



- **400lb** capacity
- Dual front castors
- New electronics with 24/7 capability
- New hand operated base opening

| PL182




| DIMENSIONS

Key dimensions of PL182	
maximum capacity	400Lb / 182Kg
base legs length	45.1"
base closed external width	25.3"
based open internal width	37.9"
mast height	48.1"
base height	4.5"
min. lifting height	25.6"
max. lifting height	69.5"
net weight	99Lb / 45Kg
shipping weight	108Lb / 49Kg


* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE AND ARE SUBJECT TO CHANGE WITHOUT NOTICE



2920 PACIFIC DR., NORCROSS, GA 30071
WWW.BESTCARELLC.COM
678-679-6690



TOP VIEW OF PL182



SIDE VIEW OF PL182

BESTLIFT™

NEW!

PL228 | FULL BODY PATIENT LIFT

The BestLift™ PL228 is the perfect blend of form and function for patient transfers. With an overall base length of 44.9", maximum width of 42.9", and base height of 4.5", it is the perfect solution for transferring bariatric patients in small rooms with low beds.



- **500lb** capacity
- Dual front castors
- New electronics with 24/7 capability
- New hand operated base opening

| PL228




| DIMENSIONS

Key dimensions of PL228	
maximum capacity	500Lb / 228Kg
base legs length	45.1"
base closed external width	25.3"
based open internal width	37.9"
mast height	48.1"
base height	4.5"
min. lifting height	25.6"
max. lifting height	69.5"
net weight	99Lb / 45Kg
shipping weight	108Lb / 49Kg


* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE AND ARE SUBJECT TO CHANGE WITHOUT NOTICE



2920 PACIFIC DR., NORCROSS, GA 30071
WWW.BESTCARELLC.COM
678-679-6690



TOP VIEW OF PL182



SIDE VIEW OF PL182

BESTLIFT™

PL400 | FULL BODY PATIENT LIFT

The BestLift™ PL400 is one of the latest ProCare products. Ergonomically designed to reduce and eliminate back injuries among professional caregivers. Providing comfort and dignity for users.




- **400lb** capacity
- Lift from floor and over high beds
- 6-point spreader bar for extra comfort
- Full electric lift with optional spare battery
- Optional Scale

| DIMENSIONS


Key dimensions of PL400	
maximum capacity	400Lb / 182Kg
external base width (closed)	26.9"
internal leg width (opened)	38.5"
upright mast height	50.6"
overall base length	47.6"
base height	4.9"
minimum lifting height with spreader bar	28.7"
maximum lifting height with spreader bar	72.2"
net weight	97lb / 44kgs
shipping weight	121lb / 54.9kgs

* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE
AND ARE SUBJECT TO CHANGE WITHOUT NOTICE

| PL400



TOP VIEW OF PL400



SIDE VIEW OF PL400

BESTLIFT™

PL500 | FULL BODY PATIENT LIFT

The BestLift™ PL500 is one of the latest ProCare products. Ergonomically designed to reduce and eliminate back injuries among professional caregivers. Providing comfort and dignity for users.




- **500lb** capacity
- Lift from floor and over high beds
- 6-point spreader bar for extra comfort
- Full electric lift with optional spare battery
- Optional scale.

| DIMENSIONS


Key dimensions of PL500	
maximum capacity	500Lb / 228 Kg
external base width (closed)	26.9"
internal leg width (opened)	38.5"
upright mast height	50.6"
overall base length	47.6"
base height	4.9"
minimum lifting height with spreader bar	28.7"
maximum lifting height with spreader bar	72.2"
net weight	98lb / 44.5kgs
shipping weight	122lb / 55.3kgs

* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE
AND ARE SUBJECT TO CHANGE WITHOUT NOTICE

| PL500



TOP VIEW OF PL500



SIDE VIEW OF PL500

BESTLIFT™

PL600 | FULL BODY PATIENT LIFT

The BestLift™ PL600 is one of the latest ProCare products. Ergonomically designed to reduce and eliminate back injuries among professional caregivers. Providing comfort and dignity for users.




- **600lb** capacity
- Lift from floor and over high beds
- 6-point spreader bar for extra comfort
- Full electric lift with optional spare battery
- Reinforced foot pedal
- Optional scale.

| DIMENSIONS


Key dimensions of PL600	
maximum capacity	600Lb / 273 Kg
external base width (closed)	26.4"
internal leg width (opened)	39.2"
upright mast height	51.3"
overall base length	47.3"
base height	4.9"
minimum lifting height with spreader bar	26.3"
maximum lifting height with spreader bar	68.7"
net weight	123 lb / 56 kgs
shipping weight	139 lb / 63 kgs

* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE
AND ARE SUBJECT TO CHANGE WITHOUT NOTICE

| PL600



TOP VIEW OF PL600



SIDE VIEW OF PL600