

Bariatric Platform Walker/Crutch Attachment

HCPCS: E0154



FIGURE A



- Designed for those who are unable to grip the walker
- Cradles forearm with soft vinyl padding for comfortable weight bearing
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker
- Adjustable strap securely holds arm in place
- Attachment for Underarm Crutch (Figure A) for item #10105 must be ordered separately (Cannot be used with Bariatric model)
- Bariatric Weight Capacity 500 lb

PRODUCTS	Description
10105HD-2	Walker Platform Attatcmt,HeavyDuty,2/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	500 lb
Height	40.5" - 55"
Carton Shipping Weight	9 lb
Weight	4 lb