

Go-Lite Bariatric Steel Rollator

HCPCS: E0149/E0156



FIGURE A



FIGURE B

- 8" casters with soft-grip tires are ideal for indoor and outdoor use
- Soft padded oversized seat
- Comes with large 20.5" (W) x 4.5" (H) x 9" (D) basket that can be mounted under seat
- Tool-free removable padded backrest for comfort
- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure A)
- Handles easily adjustable in height with self-threading knob (Figure B)

PRODUCTS

Description

10215BL-1	Rollator,Overszd w/LpBk BL,1/cs
10215RD-1	Rollator,Ovrszd,w/LpBk Red1/cs

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	18" (W) x 13" (D)
Handle Height	34" - 38.5"
Weight Capacity	500 lb
Seat To Floor Height	22"
Carton Shipping Weight	30.5 lb
Weight	26 lb