

Adjustable Height Rollator, 6" Casters

HCPCS: E0143+E0156,E0143+E0156



- Universal height adjustment on frame allows seat height to be adjusted from 18" to 22" in 1" increments
- Comes with a comfortable, durable padded seat
- Removable, hinged, padded backrest can be folded up or down
- 6" casters with soft grip tires ideal for indoor and outdoor use and make the ride smooth and stable
- Ergonomic handles are easy-to-grip and are adjustable in height with self-threading knob

PRODUCTS	Description
RTL10261BL	ADJ Height Rollator,6" Wheels,Blue,1EA
RTL10261RD	ADJ Height Rollator,6" Wheels,Red,1EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	14 lb
Material	Aluminum
Seat Dimensions	13" (W) x 12" (D)
Handle Height	29.5" - 38"
Weight Capacity	300 lb
Seat To Floor Height	18" - 22"
Casters	6"
Width Between Arms	17"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]