

Deluxe, Tilt-Top Overbed Table

HCPCS: E0274



FIGURE A



FIGURE B

- Easily assembled
- Walnut, wood grain low-pressure laminate top
- Top can be tilted 33° in either direction (Figure A)
- Top can be raised or lowered in infinite settings between 29.5"- 46"
- Swivel casters allow for easy maneuverability
- Chrome-plated steel "H" base provides security and stability
- Tabletop locks securely when height adjustment handle is released It can be raised with the slightest upward pressure
- For proper assembly, do not remove bottom screw (Figure B)

PRODUCTS	Description
13008	Over Bed Table,Tilt,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	50 lb
Base Dimensions	15.5" (D) x 26.5" (W)
Carton Shipping Weight	23 lb
Weight	21 lb