

Trapeze Base

HCPCS: E0940




FIGURE A



- For use when metal headboards are not available
- For use with home-style beds
- Attractive Brown-Vein finish (Figure A)
- Easy- to-assemble No tools necessary
- Folds easily for storage and transport
- Adjustment levers on legs

PRODUCTS

Description

13017BV	Trapeze Base, Brown Vein,1/ea
---------	-------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	250 lb
Carton Shipping Weight	27.5 lb
Weight	24 lb