

Free-Standing Silver Vein Trapeze

HCPCS:



- Designed to assist larger individuals to change positions while in bed and aid in the transfer from bed to chairs
- Easy-to-assemble (Requires only a drill, rubber mallet and socket with wrench)
- The overhead boom swivels from 0° to 180°
- Chain is adjustable from edge of headboard to the end of boom by sliding to desired position

PRODUCTS	Description
13039SV	BariatricFreeStandingTrapeze SV,1/ea
SPECIFICATIONS	Description
Weight Capacity	1000 lbs
Boom Length	36"
Carton Shipping Weight	110 lbs.
Weight	100 lbs.