

Digital Scale

HCPCS:



- Determines the weight of a person suspended from a portable patient lift mechanism
- Easy to set up and operate
- Operates with 1.9V battery
- Battery power is saved by automatic shut-off feature
- Can be used with most manufacturers lifts
- 13046-HD works with Drive's Full Body Lifts: 13240, 13242 (discontinued), 13244, 13245 (discontinued), 13023 all versions, 13010SV (discontinued)

PRODUCTS	Description
13046-HD	Digital Scale,HD 600 LB Lift,1/ea
SPECIFICATIONS	Description
Warranty	90 days
Weight Capacity	600 lbs