

Bariatric Overbed Table

HCPCS:



- In the locked position, the table can support 500 lb (Figure A)
- Structurally reinforced with a steel welded frame
- H Base Opal Frame without vanity
- Frame extends virtually the entire length of the top
- Laminate table top height adjusts from 30" to 45" in 1" increments
- 18" x 40" to reach the center line of most bariatric beds
- 2" swivel casters make for easy maneuverability
- Powder-coated base and column, with additional powder-coat finishes

PRODUCTS	Description
13080-HOBO	OBT-Bariatric-BO-H-Opal-SpclOrd
13080-HODM	OBT-Bariatric-DM-H-Opal-SpclOrd
13080-HOEO	OBT-Bariatric-EO-H-Opal-SpclOrd
13080-HOFM	OBT-Bariatric-FM-H-Opal-SpclOrd
13080-HONT	OBT-Bariatric-NT-H-Opal-SpclOrd
13080-HOQO	OBT-Bariatric-QO-H-Opal-SpclOrd
13080-HOSO	OBT-Bariatric-SO-H-Opal-SpclOrd
13080-HOWC	OBT-Bariatric-WC-H-Opal-SpclOrd

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	500 lb
Height	30" x 45"
Base Dimensions	18" (D) x 42" (W)
Carton Shipping Weight	75 lb
Weight	70 lb