

Seat Lift Chair Table

HCPCS:



FIGURE A



- Designed for use with a lift chair, standard recliner, or couch
- Can be used for many activities such as eating, writing, crafts, working on the computer, and reading
- Tabletop pivots 360° to allow for easy access and departure (Figure A)
- Right or left positioning with 180° rotation accommodates seating near any wall
- Tabletop pivots 360° at the touch of a hand
- Built-in magazine rack provides added convenience
- Teak wood grain laminate finish

PRODUCTS	Description
13085LN	Overbed Table,LiftChair,Left,1/ea
13085RN	Overbed Table,LiftChair,Right,1/ea
SPECIFICATIONS	Description
Warranty	One Year Limited
Weight Capacity	50 lb
Carton Shipping Weight	44 lb
Tabletop Adjustment	23.5" - 33"
Weight	41 lb