

Med-Aire Plus 8" Alternating Pressure and Low Air Loss Mattress System

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Active combination therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent "bottoming out" and provide up to 24 hours of power outage protection
- 8" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable 4 way stretch PU cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS	Description
14029	MedAire8 LowAirLoss.w/Pump system,1/cs
14029-84	Med Aire 8" Low Air Loss APM w/pump 84
SPECIFICATIONS	Description
Warranty	18 Month Limited
Weight Capacity	450 lb
Control Unit Airflow	8 LPM
Control Unit Alarms	Visual/Audible
Control Unit Cycle Time	10, 15, 20, 25 Minutes
Control Unit Type	Compressor
Mattress Cover Material	Fluid-Resistant Stretch
ACCESSORIES	Description
14427QF	Quick Fill Pump w/adapters, 1 EA/12 CV, EA
14427QF	Quick Fill Pump w/adapters, 1 EA/12 CV, EA