

Check and Go Oxygen Analyzer

HCPCS:



- 3-digit LCD display
- Ultra-light weight & stable
- Maintenance free
- Extended battery life
- Easy calibration
- Auto off

PRODUCTS	Description
18580	O2 analyzer & sensor,disp 1ea
SPECIFICATIONS	Description
Warranty	Limited 2 Year
Accuracy and Linearity	1% of full scale constant temperature, R.H and pressure when calibrated at full scale
Operating Temperature	60â F - 104â F (15â C - 40â C)
Resolution	0.1% Oxygen
Response Time	90% of final value < 15 seconds @ 770 F (25â C)
Warm-up Time	None
Weight	Approx 3 oz