

Aluminum Rehab Shower Commode Chair with Four Rear-locking Casters

HCPCS:



FIGURE A



- Versatile 3-in-1 design enables this model to easily be positioned over a standard toilet, used as a portable self-contained commode or shower chair.
- Simply and quickly adjust the height to achieve the ideal seated position.
- Padded flip-back armrests provide superior support and create a comfortable, non-slip surface for a safe and secure grip when sitting or standing.
- Easy to empty commode bucket to reduce splashes and spills. Use commode liners for even quicker clean up and less mess!
- 5" locking casters allow for extra security when the unit is not in motion.
- Durable and easy-to-clean padded seat (18 x 17) and backrest.
- Seat opening allows for unobstructed access to the perineal area for easy and effective cleansing.
- Swing-away detachable footrests with height and angle adjustments.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 275 lbs.
- Simple, tool-free assembly.

PRODUCTS	Description
NRS185007	Shower&Commode Chair,1/ea
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	18" (W) x 17" (D)
Weight Capacity	275 lb
Seat To Floor Height	20" - 22"
Carton Shipping Weight	40 lb
Weight	35.25 lb
Back Height	19 "