

Bariatric Anterior Safety Roller

HCPCS: E1399-1,E1399-1




FIGURE A



FIGURE B




FIGURE C



FIGURE D

- Welded steel frame
- Height adjustable in 1" increments
- Width adjustable (Figure A)
- Adjustable-brake spring tension
- Patient is situated in the center of the Safety Roller, which provides maximum support and improved posture
- Skid-resistant, non-marring, easy-rolling neoprene wheels Ideal for indoor and outdoor use
- Folds easily and stands on its own in the folded position (Figure B)

PRODUCTS	Description
CE 1000 XL	BariatricSRoller,500lb cap,Blk, 1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Carton Shipping Weight	24 lbs.135 lbs.
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
CE 1250	OxygenTank Holder,1/bx, BX
CE 1400	Dial-A-Speed Tab,1pr/bx, BX
CE 1500	Caster Lock Bracket,1pr/pkg, EA
CE 1286	Seat, Anterior Safety Rollers, 1/ea, BX
CE 1270	Extended Uprhts,Roller,1pr/bx, BX
CE 1035 FP	ForearmPlatform&MountingBracket,1pr/bx, PR