

Bariatric Offset-Handle Canes

HCPCS: E0100



FIGURE A




FIGURE B

- Manufactured with steel tubing
- Easy-to-use, dual-button height adjustment with locking ring Locking ring prevents rattling
- Comfortable vinyl hand grip with wrist strap
- Model 10305 Weight Capacity 500 lb
- Model 10318 Weight Capacity 600 lb

PRODUCTS	Description
10318-6	Offset Cane,HD,Extra Long,6/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Handle Height	29.75" - 37.75" 37" - 46"
Weight Capacity	(10305-6) 500 lb, (10318-6) 600 lb