

Cellulose Fiber Mattress

HCPCS: E0272



- 12 oz thermal-bonded polyester fiber
- Two-sided
- Inverted seam with zipper cover
- Fluid proof
- Non-allergenic
- Easy cleaning with mild soap and water

PRODUCTS	Description
15007	Mattress 80" Fiber-core, 1/EA
SPECIFICATIONS	Description
Weight Capacity	350 lbs
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]