

Clamp-On Tub Rail

HCPCS: E0246



- Durable, white powder-coated steel construction
- 14.5" - 17" height adjustment (for RTL12036-ADJ only)
- Width adjusts from 3" - 7"
- Stainless steel locking mechanism and non-wearing rubber pads provide scratch-proof security
- Fits most tubs (Tighten carefully when using with Fiberglass tubs)
- Tool-free installation

PRODUCTS	Description
12036	Tub Rail,Clmp-On,White,1/ea, RTL
RTL12036-ADJ	Tub Rail,Clmp-on,Adj.Ht.White,RTL,1 EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	300 lb
Height	14.5"- 17" 14.5"
Carton Shipping Weight	8.5 lbs 9 lb
Weight	6.5 lb 7.0 lb