

Clever-Lite Walker, Adult, with 8" Casters

HCPCS: E0143+E0156



FIGURE A



FIGURE B



FIGURE C

- Front wheels can be set in either "swivel" (Figure A) or "fixed" (Figure B) position
- Special loop lock made of internal aluminum casting operates easily and ensures safety
- Sturdy 1" diameter anodized aluminum construction
- Comes with flip-up seat with built-in carry handle
- Allows individual to be seated or use as traditional walker
- Folds with Drive-innovated side-paddle release

PRODUCTS	Description
10243	Clever-Lite Walker,w/8 wheel,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	19" (W) x 11.25" (D)
Handle Height	35.5" - 39.5"
Weight Capacity	300 lb
Carton Shipping Weight	21 lb
Weight	18 lb
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]