

Clever-Lite LS, Adult Walker

HCPCS: E0143



FIGURE A



FIGURE B

- Rear glide push-down brakes (Figure A)
- Sturdy 1" diameter anodized, extruded aluminum construction
- Comes with flip-up seat with built-in carry handle
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame
- Easily folds with patented dual-lever side-paddle release (Figure B)
- Soft, flexible backrest provides comfort and security while seated
- 5" fixed front wheels

PRODUCTS	Description
10233	Clever-Lite LS Walker, w/ 5",1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	18.75" (W) x 11" (D)
Handle Height	29.5" - 39"
Weight Capacity	300 lb
Carton Shipping Weight	15 lb
Casters	5"
Weight	11 lb
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]