

Comfort Grip Cane



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Manufactured with sturdy, extruded aluminum tubing
- Comes Standard with premium wrist strap
- Easy-to-use, one-button height adjustment with locking ring Locking ring prevents rattling
- Available in six stylish and fun patterns and colors: Anchors, Patriotic, Lemons, Rose Gold, Forest Green, and Graphite
- Handle height adjusts from 30" to 39"
- Weight Capacity 300 lb
- Approx. user height 5'-6'5"

PRODUCTS

Description

RTL10336AN	Cane T-Handle Erg, Anchor 1/ea,6/cv
RTL10336FG	Cane T-Handle Erg, Green 1/ea,6/cv
RTL10336GR	Cane T-Handle Erg, Graphite 1/ea,6/cv
RTL10336LM	Cane T-Handle Erg, Lemons 1/ea,6/cv
RTL10336MA	Cane,T-Handle Erg, Marble. 1/EA,6/CV
RTL10336MN	Cane,T-Handle Erg, Mech. Navy. 1/EA,6/CV
RTL10336PA	0ane T-Handle Erg, Patriotic 1/ea,6/cv
RTL10336RG	Cane T-Handle Erg, Rose Golds 1/ea,6/cv
RTL10336SF	Cane,T-Handle Erg, Spring. 1/EA,6/CV

SPECIFICATIONS

Description

Warranty	Limited Lifetime
----------	------------------