

Cruiser III Wheelchair

HCPCS: K0003,K0003



- Built-in seat rail extensions and extendable upholstery easily adjust seat depth from 16" to 18"
- Weighs under 36 lb (excluding front riggings)
- Carbon steel frame with silver vein finish
- Removable flip-back arms allow for easy transfer
- Frame style eliminates seat guides and allows for custom back inserts and accessories
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- Composite, Mag-style wheels are lightweight and maintenance free
- Dual axle provides easy transition of seat height to hemi-level
- 8" front casters are adjustable in three positions
- Padded armrests provide additional comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Comes with push-to-lock wheel locks

PRODUCTS	Description
K316ADDA-ELR	Cruiser III, 16", Ht. Adj. FB DDA, ELR
K316ADDA-SF	Cruiser III, 16", Ht. Adj. FB DDA, SF
K316ADFA-ELR	Cruiser III, 16", Ht. Adj. FB DFA, ELR
K316ADFA-SF	Cruiser III, 16", Ht. Adj. FB DFA, SF
K316DDA-ELR	Cruiser III, 16", FB DDA, ELR
K316DDA-SF	Cruiser III, 16", FB DDA, SF
K316DFA-ELR	Cruiser III, 16", FB DFA, ELR
K316DFA-SF	Cruiser III, 16", FB DFA, SF
K318ADDA-ELR	Cruiser III, 18", Ht. Adj. FB DDA, ELR
K318ADDA-SF	Cruiser III, 18", Ht. Adj. FB DDA, SF
K318ADFA-ELR	Cruiser III, 18", Ht. Adj. FB DFA, ELR
K318ADFA-SF	Cruiser III, 18", Ht. Adj. FB DFA, SF
K318DDA-ELR	Cruiser III, 18", FB DDA, ELR
K318DDA-SF	Cruiser III, 18", FB DDA, SF
K318DFA-ELR	Cruiser III, 18", FB DFA, ELR
K318DFA-SF	Cruiser III, 18", FB DFA, SF
K320ADDA-ELR	Cruiser III, 20", Ht. Adj. FB DDA, ELR
K320ADDA-SF	Cruiser III, 20", Ht. Adj. FB DDA, SF
K320ADFA-ELR	Cruiser III, 20", Ht. Adj. FB DFA, ELR
K320ADFA-SF	Cruiser III, 20", Ht. Adj. FB DFA, SF

PRODUCTS	Description
K320DDA-ELR	Cruiser III, 20", FB DDA, ELR
K320DDA-SF	Cruiser III, 20", FB DDA, SF
K320DFA-ELR	Cruiser III, 20", FB DFA, ELR
K320DFA-SF	Cruiser III, 20", FB DFA, SF

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	35 lb 36 lb 38 lb
Weight Capacity	300 lb 350 lb
Height	35"
Seat To Floor Height	17.5" - 19.5"
Carton Shipping Weight	43 lbs or less
Back Height	16"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES	Description
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS4S2432	Caster with 11mm Bearings, 1/ea, EA
STDS819	Anti Tippers,Cruiser,1pr/bx, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS801	Brake Handle Ext.,8", 1/ea, EA
WASR	Limb Support, Right,1/cs, CS
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS8HK-K3	Super Hemi Kit, Cruiser 3 & X4, 1/ea, EA
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS