

Deluxe 3 Wheel Aluminum Rollator, 7.5" Casters

HCPCS: E0143



- Comes standard with basket, tray and pouch (Figure A)
- Basket conveniently folds for easy storage (Figure B)
- Ships in ultra-compact carton Comes with frame folded and compact Pull lightweight, aluminum trombone-style frame to lock in place and secure with lock knob
- 7.5" casters for indoor or outdoor use
- Handles easily adjust in height with self-threading knob (Figure C)
- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure D)
- Easy, one-hand folding

PRODUCTS**Description**

10289BL	3 Whl Alum Rollator,Blue,1/ea
10289RD	3 Whl Alum Rollator,Red,1/ea

SPECIFICATIONS**Description**

Warranty	Limited Lifetime
Handle Height	31" - 38"
Weight Capacity	300 lb
Carton Shipping Weight	16 lb
Overall Width	24.7
Weight	14 lb
Width Between Arms	19.7"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]