

Deluxe Bariatric Drop-Arm Commode

HCPCS: E0168



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy-to-maintain
- Width Between Legs: 26"
- Large, durable snap-on seat
- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode
- Padded armrests provide comfort and stability
- Comes complete with 12 qt commode bucket and cover that removes from front
- Pail easily slides in from the front of the commode
- Weight Capacity 1000 lbs
- Made using recycled materials

PRODUCTS

Description

11135-1	Commode,Oversized DropArm,1/cs
---------	--------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	23.25" (W) x 18.5" (D)
Weight Capacity	1000 lbs
Carton Shipping Weight	26 lb
Height Adjustment	17.5" - 22"
Weight	21 lb