



Deluxe Folding Travel Walker with 5" Wheels

HCPCS: E0143









FIGURE A

FIGURE B

FIGURE C

FIGURE D





FIGURE E

FIGURE F

- Easily folds tool-free to ½ the size of a traditional walker (Figure A B) and can be stored in a "Ready Set Go" Pouch (Figure C)
- Legs are attached with strong elastic cords so legs can be folded (Figure D)
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces
- Two side-storage pouches with large pockets for personal items are held in place by hook and loop fastener tabs to both side frames

| PRODUCTS | Description |
|------------------------|-----------------------------------|
| RTL10263KDR | Walker w/ 5" Whls,Travel,Red,1 EA |
| SPECIFICATIONS | Description |
| Warranty | Limited Lifetime |
| Weight Capacity | 300 lb |
| Height | 29.5" - 36.5" |
| Carton Shipping Weight | 8.5 lb |
| Overall Width | 26.5" |
| Weight | 6 lb |
| Wheel Size | 5" |