

Deluxe Folding Walker, Two Button

HCPCS:



FIGURE A



- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Unique U-shaped frame design creates greater clearance
- The patient is able to bring the walker closer to assist in standing because of the improved crossframe (Figure A)
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight

PRODUCTS	Description
10200-1	Walker,2Btn.Adult,1/ea
10200-4	Deluxe Two Button Walker, 4/CS
10201-4	Deluxe Two Button Junior, 4/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	350 lb
Carton Shipping Weight	26 24 lb
Overall Width	20"
Weight	5.5 lb
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]