

Deluxe Steel Drop Arm Commode

HCPCS: E0165,E0165



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

Drop arm for ease of transfers, lightweight commode with increased width between arms.

PRODUCTS

Description

11125KD-1	Commode Drop Arm,Knock Down, 1/cs
-----------	-----------------------------------

SPECIFICATIONS

Description

Product Weight	18 lb
Seat Dimensions	13.75" - 16.5"
HCPCS	E0165
Height Adjustment	20" - 24"
Width Between Arms	19.75"
CA Prop 65	no