

Duet Rollator/Transport Chair, 8" Casters

HCPCS: E0143+E0156



- Combines the features of a rollator and transport chair in one unit An individual can ambulate independently or be pushed safely by a caregiver
- To use as a transport chair the Duet has locking flip-up armrests (Figure A - Figure C) and a flip-down footrest (Figure D) to aid the user A safety lock release is used to store the arms and a hook and loop fastener strap is used to store the footrest while using the Duet for ambulation
- The tool-free two-position contoured padded backrest is reversible depending on whether you are using the Duet as a rollator or in transport chair mode
- Comes with padded seat with carry pouch for added privacy and security of personal items Inside the pouch is a handle to aid with the folding
- The folding frame has a lock-out feature to ensure the unit stays open while in use and ensures the unit will stay folded during transport
- 8" casters are ideal for indoor and outdoor use
- Handles are adjustable in height, with self-threading knob (Figure E)

PRODUCTS	Description
795B	Duet Rollator/TransportChair,Blue,1/cs
795BK	Duet Rollator/TransportChair,Black,1/cs
795BU	Duet Rollator/TransportChair,Burgdy,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	13.5" (W) x 13" (D)
Handle Height	31.5" - 37"
Weight Capacity	300 lb
Seat To Floor Height	21"
Carton Shipping Weight	25 lb
Weight	20 lb
Back Height	35"