

Lightweight Expedition Aluminum Transport Chair

HCPCS: E1038



FIGURE A




FIGURE B

- Weighs only 19 lbs
- Available in blue or red frames
- Folds flat for easy transport (Figure A)
- Handy lever on chair back makes it easy and convenient to fold back down (Figure B)
- Comes with seat belt
- Companion activated wheel locks
- 12" rear "flat-free" wheels

PRODUCTS	Description
EXP19LTBL	Expedition X-Light Blue TransChair,1/ea
EXP19LTRD	Expedition X-Light Red TransChair,1/ea
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	19" (W) x 16" (D)
Weight Capacity	300 lbs
Seat To Floor Height	19"
Carton Shipping Weight	24 lbs
Overall Width	22" (open) 10.75" (closed)
Weight	19 lbs
Back Height	18"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]