

Folding Steel Commode

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Easily opens and folds
- Folds flat for convenient storage and transportation
- Tool-free set-up
- Durable plastic snap on seat and lid
- Comes complete with 7.5 qt commode bucket, cover and splash shield
- Grey powder coated steel welded construction increases strength and durability
- Plastic arms for added comfort
- Limited Lifetime Warranty
- Piggy back label supplied on product hangtag
- Made using recycled materials

PRODUCTS

Description

11148-1	Commode, Folding,1/cs
---------	-----------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	13.5" (W) x 15" (D)
Weight Capacity	350 lbs
Carton Shipping Weight	13.5 lbs
Weight	13.5 lbs