

Gel "E" Skin Protection 3" Gel/Foam Wheelchair Cushion

HCPCS: E2603




FIGURE A



FIGURE B

- Pressure redistributing cushion designed to assist in the prevention, treatment and management of pressure injuries
- Viscous gel bladder is encased in polyurethane foam shell and provides optimal pressure redistribution, support and comfort
- Cover top is a urethane-coated nylon that is low-shear, vapor permeable and water-resistant
- Cover base is waterproof vinyl for durability
- Carry handle included for convenient transport
- Meets flammability standard TB117

PRODUCTS	Description
14878	Gel "E" W/C Cushion, 18x18x3 ,1/ea, 4/cv
14886	Gel "E" W/C Cushion, 18x16x3, 1/cv, 4/ct
14891	Gel "E" W/C Cushion, 16x16x3, 1/bx, 4/cv
14894	Gel "E" W/C Cushion, 20x16x3, 1/cv, 4/ct
14902	Gel "E" W/C Cushion, 22x18x3, 1/cv, 4/ct
SPECIFICATIONS	Description
Warranty	18 Months Limited
Weight Capacity	Up to 20" 275 lbs 22" and up 500 lbs
Length	Fluid Resistant Nylon Top Waterproof Vinyl Bottom
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]