

Gel "E" 2" General Use Gel/Foam Wheelchair Cushion

HCPCS: E2601



FIGURE A



FIGURE B




FIGURE C

- The Medicare approved General Use 2" gel-foam wheelchair cushion will provide comfort and support for the wheelchair bound patient
- High-Density foam encapsulates a viscous gel bladder that provides proper immersion and envelopment for pressure redistribution
- The cover is fluid resistant low-shear nylon zippered and removable with a durable vinyl bottom
- Carry handle included for convenient transport
- Product Dimensions: 18" (W) X 16" (D) X 2" (H)
- Meets flammability standard TB117

PRODUCTS	Description
14888	W/C Cushion,2in Gel/Foam,18x16x2in,1/cs
SPECIFICATIONS	Description
Warranty	1 Year Limited
Weight Capacity	250 lbs
Cover Material	Fluid-Resistant Nylon Waterproof Vinyl Bottom
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]