

## Half Length Bed Rails

HCPCS: E0305




FIGURE A




FIGURE B

- Unique design prevents individuals from being lodged between bars
- Ideal for maximum protection when bed is elevated
- Fits most hospital beds with steel frame (Figure A)
- Attractive brown-vein finish is easy-to-maintain and less institutional looking (Figure B)
- Constructed of 1" steel
- Not for use with "Pan-Style" bed surface
- Easy-to-install

| PRODUCTS       | Description                             |
|----------------|-----------------------------------------|
| 15201BV        | Bed Rail, No-Gap, Half, BV, 1 pair 1/cs |
| SPECIFICATIONS | Description                             |
| Warranty       | Limited Lifetime                        |
| ACCESSORIES    | Description                             |
| 14198BV        | Clamp,15201BV No Gap Rail, 1/ea, BX     |