

Hinged Toilet Seat Riser

HCPCS: E0244-1,E0244-1



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Raises the height of the toilet 3.5" for individuals that require assistance sitting on and standing from the toilet
- Hinge design allows riser to be moved out of the way when not in use (Figure A)
- Heavy-duty molded plastic construction provides strength and durability
- Lightweight, portable, and easily attaches to toilet bowl
- Easy-to-install on most standard toilets
- Two models to accommodate either standard or elongated (Figure B) toilet bowls
- Made using recycled materials

PRODUCTS

Description

RTL12607	Standard Hinged Toilet Seat, 1 RTL/EA
RTL12608	Elongated Hinged Toilet Seat, 1 RTL/EA

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	13.5" (W) x 17.5" (D) x 3.5" (H) 14" (W) x 19.25 (D) x 3.5" (H)
Weight Capacity	250 lb
Weight	5.07 lb 5.62 lb