

Home Bed Side Helper

HCPCS:



FIGURE A




FIGURE B



FIGURE C

- Attaches to Hollywood-style metal bed frame with four legs on the floor (Figure A)
- Designed to provide assistance getting in and out of home bed
- Easy-to-install, no tools required
- Rail features four easily accessible grasping points
- Allows for normal bedding drape (Figure B)
- Adjustable leg height to fit most beds 11" wide frame for one-or two-handed use
- Durable, solid steel tubing
- Comes with 8" mounting bracket for a secure fit
- Retail Packaged

PRODUCTS

Description

15065R-1	Bed Side Helper,Retail,1/cs, RTL
----------	----------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	250 lbs
Rail Dimensions	11" (W) x 37" (H) x 11.5" (D)