

Lumbar Support General Use Wheelchair Back Cushion

HCPCS: E2611,E2611



FIGURE A



FIGURE B

- Memory-foam back cushion provides support, comfort and a tactile curve to improve posture
- Durable Masonite board prevents hammocking in sling-back wheelchair Easily removed as per patient preference
- Conforms to lumbar curve to support lower back
- Sealed liner protects foam core from moisture
- Nylon stretch cover is water resistant, non-slip vinyl back prevents the cushion from sliding and is waterproof
- Adjustable Straps easily attach to wheelchair
- Meets flammability standard TB117

PRODUCTS	Description
14889	GeneralUseBackCushion,18"x17",1/cs
14906	GeneralUseBackCushion,16inx17in,1/cs
14920	General Use Back Cushion,20x17,1/ea
8033	Premier One FoamBackCushion 22x19 ,1/ea
SPECIFICATIONS	Description
Warranty	18 Months Limited
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]