

Compartment Style Compressor Nebulizer

HCPCS: E0570-1,E0570-1




FIGURE A



FIGURE B



FIGURE C

- Includes disposable neb kit and tubing
- Attractive, Sleek Design
- Durable, Quiet System
- Piston-Driven Compressor
- Convenient Storage Compartment

PRODUCTS	Description
MQ6000	CompressorNeb,CompartmentStyle,1/ea,6/cv
SPECIFICATIONS	Description
Warranty	Limited 5 Year
Liter Flow	8 lpm
Nebulizer Medication Capacity	5ml
Noise Level	55 dBA
Operating Humidity Range	10% to 95% RH
Operating Pressure	10 psi
Particle Size	0.5 µm to 5.0 µm
Power Consumption	80 VA
Power Requirments	120V, 50-60 Hz
Storage Temperature Range	-20° to 70°C (-4° to 158° F)
Weight	4.3 lb
ACCESSORIES	Description
18254	Reusable Neb Kit 10/case, CS
MASK 001A	Adult Aerosol mask CS/50, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS