

## Quad Canes, Small Base

HCPCS: E0105



FIGURE A




FIGURE B



FIGURE C



FIGURE D

- Four point base, combined with offset handle provides additional stability and support
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting
- No assembly required
- Handle height adjusts from 30" to 39"
- Handle available in foam grip (10301-4) or vinyl (RTL1310)
- Weight Capacity 300 lb

### PRODUCTS

#### Description

|          |                                         |
|----------|-----------------------------------------|
| 10301-4  | Quad Cane,Sm.Base, Chrome,4/cs          |
| 10301F-4 | Quad Cane,Sm.Base,FoamGrp,Chrm,4/cs     |
| RTL10310 | QuadCane,Sm Base,BLK,Foam Grip,1EA,4/CS |

### SPECIFICATIONS

#### Description

|                 |                  |
|-----------------|------------------|
| Warranty        | Limited Lifetime |
| Handle Height   | 30" - 39"        |
| Weight Capacity | 300 lb           |