

Wedge Cushion with Stretch Cover

HCPCS:



FIGURE A



FIGURE B

- Wedge design helps prevent patient from sliding forward
- Firm foundation reduces pressure caused by sling-seat hammocks
- The top cover is fluid resistant, low shear stretch fabric
- The bottom cover is waterproof
- High-density foam increases stability and comfort

PRODUCTS	Description
8134	WedgeWheelchairCushion18x16x3-1 1ea,4/CS
SPECIFICATIONS	Description
Warranty	1 Year Limited
Weight Capacity	250 lbs
Cover Material	Fluid-Resistant Stretch Top Waterproof Vinyl Bottom
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]