

Portable Bath Step

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Makes getting in and out of the bath tub easier and safer
- Unique modular design allows the multiple steps to be stacked (Figure A)
- Step height is 4"
- When stacking, each additional step adds 2"
- Bath step comes with four removable covers on the surface
- When stacking, remove covers to allow steps to lock together (Figure B)
- Attractive, textured surface is slip resistant and comfortable
- Retail Packaged

PRODUCTS

Description

RTL12068	Bath Step,1 EA
----------	----------------

SPECIFICATIONS

Description

Warranty	Limited One Year
Weight Capacity	450 lb
Dimensions	19.5" (W) x 4" (H) x 16" (D)