

Aluminum Rollator, 7.5" Casters

HCPCS: E0143+E0156,E0143+E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Comes with a comfortable, durable padded seat
- Zippered pouch under seat for added privacy and security of personal items
- Removable, hinged, padded backrest can be folded up or down
- 7.5" casters with soft grip tires ideal for indoor and outdoor use and make the ride smooth and stable
- Ergonomic handles are easy-to-grip and are adjustable

PRODUCTS	Description
R728BL	Rollator, 7.5" 4Whl, Alum, Blue, EA
R728RD	Rollator, 7.5" 4Whl, Alum, Red, 1/cv
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Material	Aluminum
Seat Dimensions	14" (W) x 12" (D)
Handle Height	33" - 38"
Weight Capacity	300 lb
Seat To Floor Height	23"
Casters	7.5"
Dimensions	24 (W) x 28 (D)
Weight	14 lb
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]