

Premium Plastic, Raised, Elongated Toilet Seat with Lock

HCPCS: E0244



FIGURE A



FIGURE B

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Adds 5" to height of toilet seat
- Tool-free arms are 21" apart and can be removed or added as needed
- Wide opening in front and back for personal hygiene
- Heavy-duty molded plastic construction provides additional strength and durability
- Standard locking mechanism ensures safety
- Lightweight and portable
- No tools required for installation
- Easy-to-clean
- Fits most elongated toilets
- Made using recycled materials

PRODUCTS	Description
12013	Raised Toilet Seat w/Arms,1/cs, RTL
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	17" (D) x 15.5" (W) x 5" (H)
Weight Capacity	300 lb
Carton Shipping Weight	6.5 lb
Weight	4.5 lb
Width Between Arms	21"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]