

Premium Raised Toilet Seat with Removable Arms

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Allows individual to use existing toilet seat and lid
- Tool-free arms are 19.75" apart and can be removed or added as needed
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Lightweight and portable
- Easily attaches to toilet bowl
- No tools required for installation
- Easy to clean
- Two models accommodate MOST standard and elongated toilets
- Made using recycled materials

PRODUCTS	Description
12402	Rizer,RegularToiletsKD,w/RemArm,1/cs,RTL
12403	Rizer,ElongatedToilet,w/RemArms,1/cs,RTL
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	17" (D) x 13.75" (W) x 3.5" (H) 19.25" (D) x 13.75" (W) x 3.5" (H)
Weight Capacity	300 lb
Carton Shipping Weight	4.5 lb
Weight	3.5 lb
Width Between Arms	19.75"