

Swivel Seat Shower Stool

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- The rotating padded seat reduces twisting and reaching while in the shower
- Seat rotates 360° and locks every 90° (Figure A)
- The rotating top of the seat can be removed for easy cleaning
- Easy tool-free assembly
- Ideal for individuals who have balance or mobility issues
- Seat height can be adjusted from 16" to 22", in 1.18" increments
- Retail packaged
- Made using recycled materials

PRODUCTS	Description
RTL12061M	Swivel Seat Stool, 1 RTL, 4/CS
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	14" (W) x 14" (D)
Weight Capacity	300 lb
Carton Shipping Weight	28.65 lb
Weight	6.6 lb