

Duet Rollator/Transport Chair, 8" Casters

HCPCS: E0143+E0156



- Combines the features of a rollator and transport chair in one unit An individual can ambulate independently or be pushed safely by a caregiver
- To use as a transport chair the Duet has locking flip-up armrests (Figure A - Figure C) and a flip-down footrest (Figure D) to aid the user A safety lock release is used to store the arms and a hook and loop fastener strap is used to store the footrest while using the Duet for ambulation
- The tool-free two-position contoured padded backrest is reversible depending on whether you are using the Duet as a rollator or in transport chair mode
- Comes with padded seat with carry pouch for added privacy and security of personal items Inside the pouch is a handle to aid with the folding
- The folding frame has a lock-out feature to ensure the unit stays open while in use and ensures the unit will stay folded during transport
- 8" casters are ideal for indoor and outdoor use
- Handles are adjustable in height, with self-threading knob (Figure E)

| PRODUCTS | Description |
|----------|--|
| 795B | Duet Rollator/TransportChair,Blue,1/cs |
| 795BK | Duet Rollator/TransportChair,Black,1/cs |
| 795BU | Duet Rollator/TransportChair,Burgdy,1/cs |

| SPECIFICATIONS | Description |
|------------------------|---------------------|
| Warranty | Limited Lifetime |
| Seat Dimensions | 13.5" (W) x 13" (D) |
| Handle Height | 31.5" - 37" |
| Weight Capacity | 300 lb |
| Seat To Floor Height | 21" |
| Carton Shipping Weight | 25 lb |
| Weight | 20 lb |
| Back Height | 35" |

| PARTS | Description |
|---------------|--|
| 9501A79512 | Handlebar, Right, 795, 1/ea, EA |
| 9502M7951L | Elbow Rest, Left, 795, 2M Ser.#, 1/ea, EA |
| 9502M7951R | Elbow Rest, Right, 795, 2M Ser.#, 1/ea, EA |
| 750FORK-2M | Front Fork, 795, 2M Ser.#, 1/ea, EA |
| 9501A7953BK | FootrestAssm.Black,795 1A/2Y Ser.#,1ea, EA |
| 9501A7953R-2M | FootrestAssembly,Red,795 2M/20hSN#, 1/ea, EA |

PARTS

Description

| | |
|---------------|--|
| 9501A79511 | Handlebar, Left, 795, 1/ea, EA |
| 9501A7953B-2M | FootrestAssembly,Blue,7952M/20HSN#, 1/ea, EA |
| 9501A79510-RP | Roll Pin Fork Fastener,795,1A/2Y SN#,1ea, EA |
| 9501A79510-2M | Fork Stem Bolt, 795, 2M /20H SN#, 1/ea, EA |
| 9501A79506 | RearAxleBolt795,1A/2Y SN#,1/ea, EA |
| 9501A79506-2M | Rear Axle Bolt, 795, 2M /20H SN#, 1/ea, EA |
| 9501A79504 | Tote 795 Rollator, 1A/2Y SN#,1ea, EA |
| 750BCL-2M | HandBrakeAssy,Lt,795,2M/20H Ser.#, 1ea, EA |
| 750BRAKE-2M | Universal hand brake assy, 795, 1/ea, EA |
| 9500075005 | Folding Hinge Bolt,750/795, 4/pk, EA |
| 9501A75013 | Seat Bracket, 795 Rollator, 1/pr, EA |
| 9502M75013 | SeatBrackets&Screws,7952M/20H SN#, 1/pr, PR |
| 9501A75008 | 3 Hole Wheel,8"Rear,750N,795,1/ea, EA |
| STDS1011 | Wheel Bearing, 8mm x 22mm 1/ea, EA |