

Bariatric Steel Transport Chair

HCPCS: E1039



FIGURE A



FIGURE B

- Reinforced steel frame provides added support
- Supports individuals up to 450 lb
- Dual, reinforced steel cross braces
- 12" rear flat free wheels
- Heavy-duty, nylon reinforced upholstery with a back carry pouch
- Removable, reversible desk-length arms
- Padded armrests provide additional comfort
- Companion activated wheel locks
- Available in blue or red frames
- Swing-away footrests included

PRODUCTS

Description

BTR20-B	TransportChr.HD,20"Blue,1/cs
BTR20-R	TransportChr.HD,20"Red,1/cs
BTR22-B	TransportChr.HD,Blue,22",1/cs
BTR22-R	TransportChr.HD,Red,22",1/cs

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	20" (W) x 18" (D) 22" (W) x 18" (D)
Weight Capacity	450 lbs
Seat To Floor Height	19.5"
Carton Shipping Weight	50 lbs 54 lbs
Overall Width	24.5" (open) 9" (closed) 26.5" (open) 11.5" (closed)
Weight	45 lbs 49 lbs
Back Height	16"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

PARTS

Description

STDS1072-8	Seat Upholstery Screws, 5mm, 8/pkg, EA
STDS2K3003	HAND GRIP,BTR, 1pr, PR
STDS2K303R	Brake Assembly, BTR Right, 1/ea, EA
STDS2K303L	Brake Assembly, BTR Left, 1/ea, EA
STDS2K3002	REAR HUB CAP,BTR,6/bag, BAG

PARTS	Description
STDS2K3000	Rear Wheel for BTR22-B,12" Wheel,1/ea, EA
STDS2K3042	Front Axle Bolt & Nut for BTR Chair,1/ea, EA
STDS2K3009	8" Front Caster for BTR22, 1/ea, EA
STDS2K305R	Footrest for BTR Trnsprt Chr,Red,1/cs, CS
STDS2K305B	Footrest for BTR TransChair,Blue,1/cs, CS
STDSDSKVB	UniversalWheelChairVinylDeskArmPad,1/ea, EA
STDS2KBTR22B	BackUpholsteryForBTR22,w/Screws,1/ea, EA
STDS2KBTR22S	SeatUpholsteryForBTR22,w/Screws,1/ea, EA
STDS4S30RB	BlueArmAssyRightW/ArmPadBTR,1ea, EA
STDS4S30LB	BlueArmAssyLeftW/ArmPadBTR,1ea, EA
STDS4S30RR	RedArmAssyRight, S# 4S, BTR,1ea, EA
STDS4S30LR	RedArmAssyLeft, S#4S, BTR,1ea, EA
STDS2K3015	TIPPING SLEEVE,BTR, 1 pr, PR
STDS2K3031	REAR WHEEL BEARINGS,BTR,1/ea, EA
STDS2K3010	SLIDE TUBE, BTR,1/ea, EA