

Blue Streak Wheelchair

HCPCS: K0001-1,K0001-1



- Blue Powder-coated frame, black cross brace
- Flip-back, padded desk-length arm
- Nylon upholstery is durable, lightweight, attractive and easy to clean
- Solid Polyurethane tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Comes with calf strap
- Comes with swing-away footrests or elevating leg rests

PRODUCTS

Description

BLS16FBD-ELR	BlueStreak16"D/L FlipBackArms,Elev,1/ea
BLS16FBD-SF	BlueStreak16"D/L FlipBackArms,SwgFt,1/ea
BLS18FBD-ELR	BlueStreak18"D/L FlipBackArms,Elev,1/cs
BLS18FBD-SF	BlueStreak18"D/L FlipBackArms,SwgFt,1/cs
BLS20FBD-ELR	BlueStreak 20"D/FlipBackArms,EivLeg,1/ea
BLS20FBD-SF	BlueStreak 20"D/FlipBackArms,SwgFt,1/ea

SPECIFICATIONS

Description

Warranty	Limited 3 Year
Product Weight	41.25 lb
Weight Capacity	300 lb
Seat To Floor Height	19.5"
Carton Shipping Weight	47 lb
Back Height	16"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES

Description

STDS833N	Rear Anti Tipper w.wheels SSP&CS,1/pr, PR
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
WASL	Limb Support,Left,1/cs, CS

ACCESSORIES

Description

WASR	Limb Support, Right,1/cs, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS831	Heel Loops, 1pr, BX

PARTS

Description

STDS4S2209	Back Upholstery 16",Blue Streak,4S,1ea, EA
STDS4S2208	Back Upholstery,BS, SSP1, 4S, 1/ea, EA
STDS4S2210	Back Upholstery 20",Blue Streak,4S,1ea, EA
STDS5N2202	Saddle Assy BlueStreak,2/pack, PKG
STDS3J2416	SeatUphol.16"Cruiser3,1/ea, EA
STDS3J2420	SeatUphol.20"Cruiser3,1/ea, EA
STDS4S22X8	Cross Frame,Blue Streak 18", 4S,1/ea, EA
STDS4S2207	Side Frame,Left,Blue Streak,4S,1ea, EA
STDS5N2231	Caster Axle & Nut,BlueStreak,2 sets/pack, PKG
STDS4S22SF	Foot Rest Set,BlueStreak,4S,1set/pkg, PKG
STDS5N2201	Calf Strap BlueStreak, 1 ea/pkg, PKG
STDS5N2210	SlideTube&Bushings,BlueStreak,1/ea, EA
STDS5N22DL	L. Arm Assy. Complete BlueStreak,1/ea, EA
STDS5N22DR	R. Arm Assy. Complete BlueStreak,1/ea, EA
STDS4S2211	Left Side Guard W/Hrdwr,Blue Streak,1/ea, EA
STDS4S2212	Right Side Guard W/Hrdwr,Blue Streak,1ea, EA
STDS8M4603	BrakeTipSleeve,Cruiser,SSP, 1/ea, EA
STDS5N2219	Removable Arm Hardware, Blue Streak, 1ea, PKG
STDSDSKVB	UniversalWheelChairVinylDeskArmPad,1/ea, EA
TRHG	Hand Grips/Universal,1pr/bx, BX