

## Deluxe Fly-Weight Aluminum Transport Chair with Removable Casters

HCPCS: E1038



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Weighs only 19 lb
- Aluminum frame for lightweight strength and durability
- Easy folding and storage in closet or trunk
- Composite wheels are lightweight and maintenance free
- 8" wheels in front and rear easily remove without tools
- Deluxe back release
- Deluxe, all-aluminum rear wheel locks
- Comes with carry pocket on backrest
- Tool-free, swing-away footrests included
- Seat belt included for added safety
- Flared arms allow the back to be folded virtually flat for easy storage and transport
- Soft, plush nylon upholstery
- Convenient cup holder can easily attach to chair (sold separately Model# STDS1040S)
- Weight Capacity: 300 lb

PRODUCTS	Description
RTLFW19RW-RD	Flyweight,19in Rem Whls,Red,RTL,1 EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	19 lb
Product Dimensions	33" (L) w/ riggings x 22" (W) x 37.25" (H)
Seat Dimensions	19" (W) x 16" (D)
Weight Capacity	300 lbs
Seat To Floor Height	19"
Carton Shipping Weight	24 lbs
Folded Dimensions	33" (L) w/ riggings x 9" (W) x 37.25" (H)
Back Height	18"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]
ACCESSORIES	Description
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX